

Maximizing Efficiency in Point-of-Care Testing

i-STAT 1 Wireless expedites the availability of diagnostic information to accelerate patient care decision-making





For in vitro diagnostic use only.

Model : 300 W

Integrated Wireless Capability Optimizes the Value of Point-of-Care Testing

The benefits of *i*-STAT, seamlessly integrated with your wireless network

Easy to integrate with Electronic Medical Records (EMR)

- Test and transmit from the bedside
- Push button or auto-transmit capability
- Employ real-time sharing of test results with broad clinical teams and caregivers


Fast and accurate

- Accelerates patient care decisions
- Provides lab-quality results in minutes
- Offers comprehensive menu of tests

Flexible

- Allows clinicians to make informed decisions from bedside or anywhere EMR is available
- Empowers nurses to remain at bedside
- Enables mobile, real-time result transmission

((*i-STAT* **1** *WIRELESS*))





Streamlined Process Advances Operational Efficiency

i-STAT 1 Wireless simplifies the testing process from bedside to EMR to treatment decision



Efficient

Improves process

- Enhances predictability that test results are in the EMR in real time
- Facilitates fast, informed treatment decisions
- Reduces logistic hurdles and saves time
 - Less back-and-forth from downloader to patient bedside
 - Dock can become recharger and be placed conveniently near power jack



Ongoing Innovation in Point-of-Care Testing

i-STAT 1 Wireless* is designed for ease of use

Simple

- Fully integrated, wireless technology
- Does not require communication via Bluetooth
- Clearly displays data in large type
- Ergonomically designed, soft keys for reliable data entry

Compatible

- Designed to be compatible with 802.11b and 802.11g wireless networks
- Integrates with your existing POC data manager
- Can be interfaced with major LIS vendors

Supported

- The APOC team partners with you to ensure smooth and successful implementation
- The Technical Support Center is available 24 hours a day

To learn more about the *i-STAT 1 Wireless System*, or to learn about our other technology, process, and service innovations, contact your Abbott Point of Care Representative or visit **www.abbottpointofcare.com**

*i-STAT 1 Wireless is not available in all regions. Check with your local representative for availability in specific markets.

For *in vitro* diagnostic use only.

© Abbott Point of Care Inc. 400 College Road East, Princeton, NJ 08540 (609) 454-9000 (609) 419-9370 (fax) www.abbottpointofcare.com

i-STAT is a registered trademark of the Abbott group of companies in various jurisdictions. i-STAT Wireless Brochure 031426 Rev A 05/13



Model : 300 W



Technology | Process | Services